

## References

1. Onland W *et al.* Ethanol-lock technique for persistent bacteremia of long-term intravascular devices in pediatric patients. *Arch Pediatr Adolesc Med* 2006; 160:1049–1053.
2. Institute for Healthcare Improvement (IHI). 5 million lives campaign: How to guide: *Prevent central line infections* 2007. Accessed at [www.ihl.org/IHI/Programs/Campaign](http://www.ihl.org/IHI/Programs/Campaign)
3. Pennington CR, Pithie AD. Ethanol lock in the management of catheter occlusion. *J Parenter Enter Nutr* 1987; 11: 507–508
4. Johnston DA *et al.* Ethanol flush for the prevention occlusion of catheter. *Clinical Nutrition* 2002;11: 97-100
5. Centers for Disease Control and Prevention. Guidelines for the prevention of intravascular catheter-related infections. *MMWR* 2002; 51(RR-10):[inclusive page numbers].
6. Ball PA *et al.* Ethanol locking as a possible treatment for microbial contamination of long-term central venous catheters. *Nutrition* 2003; 19: 570
7. Dannenberg C *et al.* Ethanol-lock technique in the treatment of bloodstream infections in pediatric oncology patients with broviac catheter. *J Pediatr Hematol Oncol* 2003; 5: 616-621
8. Broom J *et al.* Ethanol lock therapy to treat central venous catheter-associated blood stream infections: results from a prospective trial. *Scand J Infect Dis* 2008; 40(5): 399-406
9. Ackoundou-N'guessan C *et al.* Ethanol lock solution as an adjunct treatment for preventing recurrent catheter-related sepsis—first case report in dialysis setting. *Nephrol Dial Transplant* 2006; 21: 3339
10. Metcalf S *et al.* Use of ethanol locks to prevent recurrent central line sepsis. *Journal of Infection* 2004; 49: 20–22
11. Opilla MT, Kirby DF, Edmond MB. Use of ethanol lock therapy to reduce the incidence of catheter-related blood-stream infections in home parenteral nutrition patients. *JPEN* 2007; 31(4): 302-305
12. Lairda, J *et al.* Complications of the ethanol-lock technique in the treatment of central venous catheter sepsis *Journal of Infection* 2005; 51: 338–341
13. Cober MP, Johnson CE. Stability of 70% alcohol solutions in polypropylene syringes for use in ethanol-lock therapy. *Am J Health Syst Pharm* 2007; 64 (23): 2480-2
14. Pomplun M *et al.* Stability of a heparin-free 50% ethanol lock solution for central venous catheters. *J Oncol Pharm Pract* 2007;13 (1): 33-7
15. Crnich M *et al.* The effects of prolonged ethanol exposure on the mechanical properties of polyurethane and silicone catheters used for intravascular access. *Infect Control Hosp Epidemiol* 2005; 26: 708-714
16. Guenu S *et al.* Mass spectrometry and scanning electron microscopy study of silicone tunneled dialysis catheter integrity after an exposure of 15 days to 60% ethanol solution. *Rapid Commun Mass Spectrom* 2007; 21(2): 229-36
17. Bell AL. *et al.* Ethanol/trisodium citrate for hemodialysis catheter lock. *Clinical Nephrology* 2004; 62 (5): 369-73
18. Takla TA, Zelenitsky SA, Vercaigne LM Effect of ethanol/trisodium citrate lock on microorganisms causing hemodialysis catheter-related infections. *Scand J Infect Dis.* 2008; 40 (5): 399-406
19. Raad I *et al.* Optimal antimicrobial catheter lock solution, using different combinations of minocycline, EDTA, and 25-percent ethanol, rapidly eradicates organisms embedded in biofilm. *Antimicrobial Agents and Chemotherapy* Jan. 2007: 78–83
20. Aslam S. Effect of antibacterials on biofilms. *Am J Infect Control.* 2008 Dec; 36 (10): S175.e9-11.
21. Sanders J *et al.* A prospective double-blind randomized trial comparing intraluminal ethanol with heparinized saline for the prevention of catheter-associated bloodstream infection in immunosuppressed haematology patients. *The Journal of Antimicrobial Chemotherapy.* 2008 62(4):809-815.
22. Balestrino D *et al.* Eradication of microorganisms embedded in biofilm by an ethanol-based catheter lock solution. *Nephrol Dial Transplant.* 2009 Apr 21. [Epub ahead of print]